

TRAX

Compression Screw System

Product Information

Breakaway Screw



Headless Screw



Headed Screw



The **TRAX Compression Screw System** offers a wide range of bone screws for osteotomy, fixation, fusion, and fracture procedures. With several screw and instrument options arranged into three kits, selecting the perfect implant for your patient is straightforward.

TS
TRAX SURGICAL


The TRAX Compression Screw System also includes a set of non-cannulated titanium screws in the small kit. These headed, break-away screws are designed for small bones/fragments, or where small screws are desired for optimum bone reconstruction, repair, osteotomies, and fusions.

**2.0mm Diameter
Break-Away Screws**

Length	Qty/Kit
11mm	4
12mm	4
13mm	4
14mm	4

Threads designed for rapid insertion into cancellous bone, and maintain a firm hold.

Pin breaks away below the screw head to minimize prominence. Lessens risk of soft tissue irritation.



Trax Surgical Inc.

75 Mill Street
Stoughton, MA 02072
Tel: 781-436-4350
www.TraxSurgical.com

P/N 68405 Rev. A 09/2023

